

***HyperShear™* HPLC and UHPLC Mixers**




***HyperShear™* HPLC and UHPLC Mixer Features:**

- Micro Flow, Low Flow, Analytical Flow, High Flow and Combo Series**
- Reduce baseline noise, increase sensitivity and improve gradient accuracy**
- UHPLC mixers are available in volumes ranging from 0.5 µL to 1.5 mL and pressure to 15,000 PSI**
- Available for In-line, Binary and Ternary formats in Stainless Steel (SS) or PEEK**
- Ideal for microbore HPLC, UHPLC and LC/MS**
- Compact design is easily integrated into any HPLC and UHPLC system**
- Increased reaction efficiency for post column derivatization**
- Decreased mixing and delay volume without sacrificing mixing efficiency**
- Mixing volume optimization is easy with interchangeable mixer cartridges**

HyperShear™ HPLC and UHPLC Mixers

ASI manufactures a wide range of Static Mixers to solve the most demanding high pressure mixing problems. ASI HyperShear Mixers incorporate a highly efficient cross-flow shearing mechanism which produces vortex shear mixing over a wide flow range. This mixing technology typically delivers between 25% to 200% better mixing efficiency compared to conventional packed bed or tortuous path mixers. HyperShear Mixers are available in a variety of housing formats including: inline, binary, and ternary inlets. Mixers are constructed from stainless steel or Biocompatible PEEK with volumes ranging from 0.5 µL to 1.5 mL.

Within a given mixer flow series, mixer cartridges are interchangeable. The ability to swap mixer cartridges within a given mixer series makes selecting the optimum mixer volume easy and economical. This is not the case with UHPLC mixers. Since we warranty the pressure rating of 15,000 PSI on the factory tested mixer assembly (housing and cartridge), UHPLC mixers are only sold as complete mixer assemblies.



Specifications

	Standard HPLC	UHPLC	Biocompatible PEEK
Pressure Rating	6,000 PSI	15,000 PSI	3,000 PSI
Female Port Geometry	10/32 Parker(1/16 CPI)	10/32 Parker(1/16 CPI)	10/32 Parker(1/16 CPI)
Wetted Materials	Stainless Steel and PEEK	Stainless Steel and PEEK	PEEK
Micro Flow Series 0.5, 1, 2, 5, 10 and 25 µL	available	available	available
Low Flow Series 50, 150 and 250 µL	available	available	available
Analytical Flow Series 350 and 500 µL	available	available	available (In-Line)
High Flow Series 800 µL, 1.0 and 1.5 mL	available	available	available (In-Line)
Combo Series 1 to 100 µL	available (In-Line)	available (In-Line)	available (In-Line)

Static Mixer Application Notes

Selecting the Right Size Mixer Cartridge

ASI offers static mixers with volumes that range from 0.5 microliters to 1,500 microliters. Choosing the right size mixer is a trade off between delay volume, mixing noise, gradient fidelity and chromatographic performance. **Please refer to the Charts on pages 64 and 65 for more information on gradient accuracy as a function of flow rate and mixing volume.**

The following observations will provide some guidelines to help choose the right size mixer.

- For any given flow rate, the more the mixing volume the better the mixing, and the lower the baseline noise.
- The smaller the mixing volume, the better the definition and sharpness of linear gradients.
- Multi-pump high pressure gradient systems typically require far less mixing volume than low pressure single pump gradient systems when running linear gradients.
- An ASI 150 µL in-line static mixer can be added in addition to the standard onboard mixer to further reduce mixing noise.

Multi-pump High Pressure Gradient Systems

Linear Gradients

If a larger mixing volume can be tolerated for a particular flow rate, the larger the volume will lower the mixing noise. The upper limits to mixing volume will be the maximum delay time that can be tolerated, and possible distortion (tailing) of the gradient at the beginning and end of the gradient. The lower limit will be defined by the amount of mixing noise that can be tolerated.

Please refer to a table, page 59 for specific recommendations.

Binary or Ternary – Steady State Composition

Always select the largest volume that will still provide an acceptable delay volume. In general, the more mixing volume, the better the mixing will be. For most pump systems, a 150 µL cartridge will provide adequate mixing.

Examples of this type of pump system include:

Shimadzu LC-10AD and LC-10ADvp
Beckman System Gold®
Gilson Model 305
Agilent Model 1100

Single-pump Low Pressure Mixing Gradient Systems

Linear Gradients

These systems generally require more mixing volume to perform linear gradients than multi-pump high pressure systems. The following will explain why this is the case. In a low pressure system the composition can only be changed once every pump stroke. Since the pump stroke volume of most pumps is 100 μL , and it takes a mixer volume that is about 3 times the batch volume to provide adequate mixing, we need 300 μL of mixer volume, at least, to do adequate mixing. More insoluble combinations may require even more mixing volume.

In general, choose the largest size mixer cartridge that will still provide an acceptable delay volume. For most applications this will be at least 350 μL .

Binary or Ternary – Steady State Composition

Always select the largest volume that will still provide an acceptable delay time. In general, the more mixing volume, the better the mixing will be. For most applications this will be at least 350 μL .

Examples of this type of pump system include:

Agilent 1090
 Perkin Elmer series 200
 TSP Spectra Vision®
 Waters model 626
 Varian Star® 9000

Mixer Cartridge Selection Guide for High Pressure Mixing

Linear Gradients - High Pressure Mixing

Larger mixing volume can be tolerated for a particular flow rate, with the larger the volume the lower the mixing noise. The upper limits to mixing volume will be the maximum delay time that can be tolerated, and possible distortion (tailing) of the gradient at the beginning and end of the gradient. The lower limit will be defined by the amount of mixing noise that can be tolerated. The following cartridge volumes are a compromise between these two limits.


Table 1.

Flow	Mixer Cartridge Volume
0-5 micro liter/min.	5 micro liter
5-10	10
10-20	25
20-150	50
150-500	150
500+	250

Static Mixers

Binary Static Mixer Comparison


Less Noise - Lower Detection Limits



HPLC System Configuration: ASI/Model 500G Gradient System ABI Model 785 UV/VIS @254 on-column Varian Star Data Acquisition	HPLC Conditions: MP: A = H ₂ O, B = H ₂ O/Acetone Flow Rate: 50 μL/min. 50% B Pressure: 1,700 PSI with ASI resistor cartridge
--	---

Figure 21


Superior Mixing with Less Gradient Dispersion



HPLC System Configuration: Above	HPLC Conditions: Above Isocratic Flow Rate: 50 μL/min. 100% A One inlet port plugged Injection: 600 nL Acetone
---	--


Figure 22

Agilent Static Mixer Optimization



<p>HPLC System Configuration: Agilent® Binary 1100 Pump Packed bed mixer (approximate volume: 500 µL) replaced with union and/or ASI HyperShear In-Line Static Mixer Varian Star Data Acquisition</p>	<p>HPLC Conditions: MP: A = H₂O, B = ACN doped with 0.01% Acetone Flow Rate: 250 µL/min. Pressure: 2,000 PSI with ASI resistor cartridge Gradient: 0.0 = off%B 7.5 = 20%B 2.5 = 10%B 12.4 = 20%B 7.4 = 10%B 12.5 = 0.0%B and hold</p>
---	---

Figure 23



<p>HPLC System Configuration: Above Agilent® Binary 1100 Pump Packed bed mixer (approximate volume: 500 µL) replaced with union and/or ASI HyperShear In-Line Static Mixer Varian Star Data Acquisition</p>	<p>HPLC Conditions: MP: A = H₂O, B = IPA doped with 0.01% Acetone Flow Rate: 1.0 mL/min. Pressure: 2,000 PSI with ASI resistor cartridge Gradient: Above</p>
---	--

Figure 24

Static Mixers

Shimadzu Static Mixer Optimization

Micro Flow Step Gradients





Figure 25

Micro Flow Step Gradients



HPLC System Configuration:

Shimadzu: LC-10ADvp
Mixer: Connected at Shimadzu Tee
Detector: Knauer 2501 @254


HPLC Conditions:

MP: A = H₂O, B = H₂O + 35% 2-Propanol + 0.003% Acetone
Flow Rate: 250 µL/min.
Pressure: 2,000 PSI with ASI resistor cartridge

Figure 26


Shimadzu Static Mixer Optimization

Analytical Flow Step Gradients



HPLC System Configuration: Shimadzu: LC-10ADvp Mixer: Connected at Shimadzu Tee Detector: Knauer 2501 @254	HPLC Conditions: MP: A = H ₂ O, B = H ₂ O + 35% 2-Propanol + 0.003% Acetone Flow Rate: 1.0 mL/min. Pressure: 2,000 PSI with ASI resistor cartridge
--	--



Figure 27



HPLC System Configuration: Shimadzu: LC-10ADvp Mixer: Connected at Shimadzu Tee Detector: Knauer 2501 @254	HPLC Conditions: MP: A = H ₂ O, B = H ₂ O + 35% 2-Propanol + 0.003% Acetone Flow Rate: 1.0 mL/min. Pressure: 2,000 PSI with ASI resistor cartridge
--	--

Figure 28

Effect of Mixer Volume on Gradient Accuracy with a Constant Flow Rate




HPLC System configuration:
 Shimadzu LC-10ADVP pumps with micro-flow modification
 Shimadzu SPD-10AVP UV-VIS detector
 Shimadzu Class-VP v5.03 software



HPLC conditions:
 Gradient: 10-90% B in 10 min. steps
 A = H₂O
 B = 0.3% Acetone in H₂O

Figure 29


Effect of Flow Rate on Gradient Accuracy with a Constant Mixing Volume



ASI 25 µL mixer, 100 µL/min.



ASI 25 µL mixer, 50 µL/min.



HPLC conditions:
 Gradient: 10-90% B in 10 min. steps
 A = H₂O
 B = 0.3% Acetone in H₂O

HPLC System configuration:
 Shimadzu LC-10ADVP pumps with micro-flow modification
 Shimadzu SPD-10AVP UV-VIS detector
 Shimadzu Class-VP v5.03 software

Figure 30